



# DUNAMIS DENTAL

One Handpiece per a Patient



## ■ COMPANY PROFILE



# DUNAMIS DENTAL

## One Handpiece per a Patient

Since its foundation in 2001, Dunamis has been dedicated to developing quality highspeed dental handpieces through in-house R&D and global partnership. With production capacity of 12,000 units per year, Dunamis is ready to support global distribution network. Dunamis will continue to grow to become a leading supplier of world's dental handpiece market.

KT15 HANDPIECE



KT15 CARTRIDGE



HANDPIECE LUBRICANT  
- KARAM SPRAY



## ■ Highspeed Handpiece



KT15-SQD



KT15-SK




KT15-SN



### ■ KT15-SQD Specifications

<b>Speed</b>	380,000 ~ 400,000 rpm
<b>Torque</b>	12gf·cm Min (STD)
<b>Chuck Type</b>	Push Button Type
<b>Chuck Force</b>	2.5 Kgf Min
<b>Water Spray</b>	4-Hole Spray
<b>Bearing</b>	Ceramic Ball Bearing
<b>Autoclave</b>	135°C Autoclave

### ■ Coupling end Configuration



Midwest 4-hole




Borden 2-hole


## ■ Highspeed Handpiece



KT15-SLK



KT15-SLN



KT15-ECN S M4



KT15-ECN S B2



# KT15-45 Surgical Handpiece



## KT15-45SQD

- Prevention of Emphysema
- Easy Third Molar treatment at 45° angle head
- Buried water pipe Ceramic ball bearing (400,000 rpm)
- Compact capsule cartridge type
- Single Water Injection
- Autoclaveable up to 135°C
- Swivel & Quick connecting coupling

### 1. Apicoectomy case      2. Odontectomy of Angulated, Impacted Third Molar



After enucleation of cyst, cut root apex portion.


Easily accessible without occurrence of emphysema.



A) After exposure of crown portion, carry out odontectomy.



B) After removal of crown portion, give indentation grooves to bone and roots for removal of remaining roots.



C) Pieces of removed third molar

## ▪ Lowspeed Handpiece



KARAM-ANG



KARAM-AM



KARAM-STR

## ▪ Highspeed Handpiece(Laboratory)



KARAM-DELIX



KARAM-LAB

## ▪ Handpiece Lubricant - Karam Spray



## Karam Spray

- Applicable to all dental handpieces
- Content Volume : 550ml
- Original Formula by Klüber Company (Germany)

[www.dunamis.co.kr](http://www.dunamis.co.kr)



904, "A"-Dong, Gwangmyeong Technopark, Haan-Ro 60, Gwangmyeong-Si,  
Gyeonggi-Do, Republic of Korea, 14322

**TEL** +82 (0)2 778-4565/6 **FAX** +82(0)2 778-4501 **E-mail** [yskim@dunamis.co.kr](mailto:yskim@dunamis.co.kr)

